

DISPOSABLE POLYETHYLENE GOWN - LARGE



Gown Technical Data:

- Made in the USA (no reliance on foreign SMS materials and gown manufacturers) resin, films, manufacturing and converting
- Single-layer polyethylene for exceptional protection
- Lightweight, non-linting
- Material: 1.5 MIL, LLDPE, sky blue
- Ideal for barrier: poly provides superior protection against food waste, liquids, or debris
- Apron-style for quick donning and pull-away removal
- Thumb-hook cuffs for a secure fit and complete arm coverage when worn beneath gloves
- Waist ties in the back for a custom fit
- One size (large) fits all
- Bulk-packaged, poly-lined 100/box (10 bags of 10 gowns)
- Minimum order: one box



Image rendered, not actual product

Vonco Overview:

- Location: Trevor, WI, USA
- Year established: 1955
- New state-of-the-art manufacturing facility built 2017
- Square footage: 88,000 (expansion planned)
- Unique capabilities: Manufacturing of liquid tight fluid bags, infection prevention, and biohazard transport bags
- ISO 9001 and ISO 13485:2016 Certified
- FDA registered facility with Class II medical device capability



Product Master Drawing

