

## POLY GOWN: LEVEL III SURGICAL - LARGE

### SPEC FOR AAMI LEVEL III SURGICAL POLY GOWN (NON-SMS)

Level III disposable polyethylene surgical gown meets liquid barrier protection at level III or higher, consistent with ANSI/AAMI PB70 for critical zone areas. Meets the Class I or Class II flammability standard per 16 CFR Part 1610.

#### Vonco Overview:

- Location: Trevor, WI, USA
- Year established: 1955
- New state-of-the-art manufacturing facility built 2017
- Square footage: 88,000 (expansion planned)
- Unique capabilities: Contract manufacturing of Medical Devices and Fluid Bags, Infection Prevention, and Biohazard Transport Bags
- ISO 9001 and ISO 13485:2016 Certified
- FDA registered facility with Class II medical device capability

#### Baseline Technical Data:

- Made in the USA (no reliance on foreign SMS materials and gown manufacturers) resin, films, manufacturing and converting
- Single-layer polyethylene for exceptional protection
- Lightweight, non-linting
- Material: 1.5 MIL, LLDPE, sky blue
- Apron-style for quick donning and pull-away removal
- Thumb-hook cuffs for a secure fit and complete arm coverage when worn beneath gloves
- Waist ties in the back for a custom fit
- One size (large) fits all
- Ideal for barrier, poly provides superior protection against liquid-borne pathogens, liquids, waste or debris
- Manufactured to meet FDA recognized American National Standards Institute/Association of Advancement of Medical Instrumentation (ANSI/AAMI)
  - Level III Moderate risk: for use during arterial blood draw, inserting an IV line, in the ER, or for trauma cases
- Minimum order required: 10,000 gowns
- Level III gowns are custom made to order – additional lead times and certifications apply
- Contact Vonco for pricing and delivery

