

AAMI Level 2 Non-Sterile Isolation Gown

Vonco Products Level 2 Non-Sterile Isolation Gowns meet the requirements of AAMI PB70:2012 for critical zone areas, and the flammability requirements of 16 CFR Part 1611.

Vonco Overview:

- Location: Trevor, WI, USA
- Year established: 1955
- New state-of-the-art manufacturing facility built 2017
- Square footage: 88,000 (expansion planned)
- Contract manufacturer of medical devices and fluid bags, infection prevention, and biohazard transport bags
- ISO 13485:2016 Certified
- FDA registered facility with Class II medical device capability

Baseline Technical Data:

- Made in the USA (no reliance on foreign SMS materials and gown manufacturers) resin, films, manufacturing and converting
- Single-layer polyethylene for exceptional protection
- Lightweight, non-linting
- Material: 1.5 MIL, LLDPE, sky blue
- Apron-style for quick donning and pull-away removal
- Thumb-hook cuffs for a secure fit and complete arm coverage when worn beneath gloves
- Waist ties in the back for a custom fit
- One size (large) fits all
- Level 2 non-sterile isolation gowns are intended to protect the wearer from the transfer of microorganisms and bodily fluids in low- or minimal-risk patient isolation situations
- Not recommended for use in a surgical setting or where significant exposure to liquid, bodily or other hazardous fluids may be expected, use in a clinical setting where Level 3 or 4 protection is warranted, or use in the presence of high intensity heat source or flammable gas
- Contact Vonco for samples, pricing and delivery



10826 250th Avenue
Trevor, Wisconsin 53179
800-323-9077
sales@vonco.com

