

EnteraLoc™ FLOW

FDA Over-the-Counter Cleared

This groundbreaking enteral device sets the bar for future direct connect delivery of nutrition — blended and non-blended



2022 FPA HIGH ACHIEVEMENT
AWARD RECIPIENT



EnteraLoc Flow is an enteral feeding connection system featuring a spouted pouch and ENFit[®] connector that provides safer enteral nutrition with superior nutritional intake and improved patient outcomes.

Turnkey Contract Manufacturing

Services include product design, prototyping, tooling and assembly line production, acquisition of materials, testing, product certifications, and complete supply chain management through filling.

Spout Insertion

Vonco produces spouted pouches from start to finish. Our equipment can easily insert interchangeable spouts into flexible pouches which helps us achieve fast turnaround.

Custom Spouted Pouch

Spouted pouches are becoming more prevalent in healthcare because of the increasingly wide variety of convenient spouts and fitments.

Lab Line Filling

Our filling partnerships give us the flexibility to test small batches to meet your specifications. Or, Vonco can drop lab filling equipment at your facility.

Direct Connect ENFit Spout


EnteraLoc Flow comes ready-to-use out of the box. Simply select the preferred meal and attach the direct connect device to the feeding tube. That's it.

Commercial Filling

Using commercial filling processes give us the flexibility to package retort foods in a format that meets your specifications. From small batches to large production runs, our linear fillers can automatically dose, fill and cap pouches.



Rapid Commercialization Process



Vonco reduces the complexity and cost of sourcing, manufacturing and distributing enteral feeding solutions by offering customized turnkey manufacturing services.



Visit www.vonco.com to learn more

Vonco Products, LLC.

10826 250th Avenue, Trevor, Wisconsin 53179

800.323.9077

sales@vonco.com

© 2023 Vonco Products, LLC. All rights reserved.